

## Gilson 832 Controller

---

[Read more on the listing...](#)

\*Controller only - No Oven\*

For controlled heating and cooling of samples from 4°C to 40°C.

---

Independent temperature control of up to two racks.

Temperature Range: 4 &degC to 40 &degC.

Increment: 1 &degC.

Accuracy:  $\pm 1.5$  &degC.

Stability:  $\pm 0.5$  &degC.

Unit of measurement: &degC, &degF, or K.

User Interface: 2 x 24-character LCD display and soft touch keypad.

Electrical Interfaces: Control cable, Gilson Serial Input Output Channel (GSIOC) ID#46.

Dimensions: 290mm x 160mm x 140mm.

Weight: 3.8 kgs.

Operating voltage: 230V.

This item is in good condition - -please contact us for more details.

\*\*Control unit only - no oven supplied.