

Agilent 1100 G1313A ALS Autosampler

[Read more on the listing...](#)

Year:2002

Variable Volume Autosampler with a sample capacity of 100 samples. Variable volume injections in the 1 to 100 μ L range are programmable.

Specification

Carryover: Typically less than 0.1 %, less than 0.05 % with external needle cleaning

Sample capacity: 100 x 2-ml vials in 1 tray, 40 x 2-ml vials in 1/2 tray, 15 x 6-ml vials in 1/2 tray (Agilent vials only)

Injection cycle time: Typically 50 seconds depending on draw speed and injection volume

Operating Voltage: 230 Volts AC.

The autosampler can be integrated into almost all HPLC systems to provide automated analysis and comes completely reconditioned in excellent, working condition - please contact us for any further information.