

## Agilent 1100 G1314A UV Detector

---

[Read more on the listing...](#)

Year 1999

UV Detector with wavelength range of 190nm - 600nm.

---

**Detection type:** Double-beam photometer

**Light source:** Deuterium lamp

**Wavelength accuracy:**  $\pm 1$  nm (Self-calibration with deuterium lines, verification with holmium oxide filter)

**Band width:** 6.5 nm typical

**Analogue outputs:** Recorder/integrator: 100 mV or, 1 V, output range 0.001 - 2 AU, one output

**Operating Voltage:** 230 Volts AC.

This item is in good condition - please contact us for any further information.