

Agilent 1100 G1367A Well Plate Autosampler

[Read more on the listing...](#)

Year:2005

Well plate autosampler supplied with 96 x 0.5ml polypropylene tray.

Features

Increased sample injection speed for high sample throughput.

Overlapped injections for increased productivity.

Minimal delay volumes for rapid gradients and fast equilibration when bypassing the autosampler after sample injection.

Flexible and convenient sample handling with different types of sample containers. 384-well plates allows you to process up to 768 samples unattended.

Specification

Injection range: 0.1 - 100 μ L in 0.1 μ L increments.

Operating Voltage: 230 Volts AC.

Can be integrated into almost all HPLC systems to provide automated analysis - various tray options available - please enquire.

Comes fully serviced with OQ - please contact us for any further information.