

Waters M 590 Isocratic Pump

[Read more on the listing...](#)

This reconditioned Waters M590 Programmable Isocratic pump is in a superb condition and has been tested and OQd.

This unit offers control of external system accessories, a wide flow rate range and is programmable to allow unattended operation. The 590 interfaces with Waters controllers for multipump gradient capabilities.

Specification

Flow rate range: 0.001 to 20.00mL/minute in 0.001mL increments

Flow readout range: 0.000 to 10.00mL/minute

Operating voltage: 230V

Please enquire for full details.